## Speichel SARS-CoV-2 (2019-nCoV)

### Antigen Test Kit (Nanocarbon Assay) v1.6

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
Ningbo Beautiful Life Medical Biotechnology Development Co., Ltd.  
Address: Room 305, Building 4, NO.7 Chuangye Avenue,  
Free Trade West Zone, Ningbo, China

## Saliva SARS-CoV-2 (2019-nCoV)

### Antigen Test Kit (Nanocarbon Assay) v1.6

#### For professional in-vitro diagnostics Only

#### PRODUCT NAME

Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay)

#### CATALOG NUMBER

NBMLRS-XG-S2

#### PACKAGE SPECIFICATION

1 Test/Box; 20 Tests/Box

#### INTENDED USE

Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay) is a lateral flow, one-step immunoassay for the qualitative detection of N antigen in Saliva sample. This product is used to obtain a visual, qualitative results and is intended for professional use.

This assay provides only a preliminary analytical test result. A more specific alternate method must be used in order to obtain a confirmed analytical result. The RT-PCR Test is the preferred confirmatory method.

#### DETECTION PRINCIPLES

Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay) is a rapid chromatographic immunoassay test device based on the principle of sandwich antibody and antigen binding.

During testing, SARS-CoV-2 N antigen, if the virus antigens present in the saliva specimen, it will bind to the antibody conjugates to form a complex. This complex will migrate upward by capillary and then be captured by the immobilized antibody coated at the Test line region (T line) on the NC membrane.

A visible black line will show up in the Test line region (T line).

The black line will not form in the Test line region (T line) if there is no virus antigen in the sample.

A virus antigen-positive saliva specimen will generate a colored line in the Test line region (T line) because of antigen binding, while a virus antigen-negative saliva specimen or a specimen containing a very low concentration will not generate a colored line in the Test line region (T line).

To serve as a procedural control, a colored line will always appear at the Control line region (C line), indicating that proper volume of specimen has been added and membrane wicking has occurred.

#### MAJOR COMPONENTS

Product Name	Major Biochemical Composition
Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay)	It consists of strips, test cassette, desiccants and aluminum foil bags. Strips consist of plastic substrates, colloidal Nanocarbon pads, cellulose films, loading pads, and absorbent pads. Colloidal Nanocarbon pads are pre-coated with Nanocarbon-labeled anti-N antigen monoclonal antibodies, and the detection line of nitrocellulose. The QC line is goat anti-mouse IgG antibody.


#### STORAGE CONDITIONS AND EXPIRATION DATE

The product should be stored at 2°C-30°C in a dry place and avoid light; it expires in 12 months from the date of its production and must not be frozen.

The product should be used in an hour after the aluminum foil bag is torn; if the aluminum foil bag is torn in an environment with high humidity, it should be immediately used. Do not use beyond the expiration date.


#### TESTING METHODS

##### Method A:




1. Collect enough fresh saliva in a single use disposable cup.
2. Open the aluminum foil bag, take out the test cassette, un-plug the cap.
3. Place the absorbing tip into the saliva specimens, let the absorbing tip to immerse in and absorb saliva fully.
4. Keep the test cassette upright and let the saliva fluids to move upward until reaching over Line C, then plug the cap back.
5. After 15 minutes, interpret the test result.

##### Method B:

1. Open the aluminum foil bag, take out the test cassette, un-plug the cap.
2. Place the absorbing tip beneath tongue for 2 minutes.
3. Keep the test cassette upright and let the saliva fluids to move upward until reaching over Line C, then plug the cap back.
4. After 15 minutes, interpret the test result.

#### EXPLANATION OF TEST RESULTS

	Positive (Infected)	One black line should be in the Control region (C), and another apparent black line should be in the Test region (T). This positive result indicates that there are SARS-CoV-2 N antigen presenting in the saliva specimen.
	Negative (Not Infected)	One black line appears in the Control region (C). No black line appears in the Test region (T). This negative result indicates that there are no SARS-CoV-2 N antigen presenting in the saliva specimen or the antigen concentration is very low.
	Invalid	Control line fails to appear. Insufficient specimen volume, incorrect procedure or technical error are the most likely reasons for Control line failure. Please repeat the test procedure. If the problem persists, please contact your local distributor.

### LIMITATIONS OF THE TEST METHOD

1. The test device provides only a preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. The RT-PCR test is the preferred confirmatory method.
2. The device is a screening assay and cannot determine either the concentration of the antigen in the saliva or the level of the virus.
3. It is possible that technical or procedural errors, as well as other interfering substances in the saliva specimen may cause erroneous results.
4. Adulterant, such as bleach and/or alum, in saliva specimen may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another saliva specimen.
5. A Positive result indicates presence of the antigen but does not indicate the amount of the virus in saliva specimen.
6. A Negative result may not necessarily indicate virus-free saliva. False Negative results may be obtained when virus is present but below the cut-off level of the test.

### NOTES

1. For in vitro diagnostic use only. For professional use only.
2. Do not use after the expiration date.
3. The test cassette should remain in the sealed aluminum foil bag and until use.
4. Saliva specimen may be potentially infectious. Proper handling and disposal measures should be adopted.
5. The used test cassette should be discarded according to Federal / State Laws and Local Regulations.

### CLINICAL PERFORMANCE

**Purpose:** This Clinical Performance Study is to determine the performance of the product compared with the competitors to obtain the clinical research analysis.

#### Method:

**Antigen Test:**

**Device:** Saliva SARS-CoV-2 (2019-nCoV) Antigen Test Kit (Nanocarbon Assay)

**Sample:** Saliva specimens were collected in sterile urine cups without addition of any preservatives from individual tester.

**RT-PCR Test:**

**Device:** EUA-Authorized ThermoFisher Scientific TaqPath RT-PCR COVID-19 combo kit + Applied Biosystems™ QuantStudio™ 12K Flex Real-Time PCR System instrument.

**Sample:** Nasopharyngeal swab specimens were collected in 3 mL viral transport media from individual tester.

The Clinical Performance Study enrolled 244 subjects to randomize the blinded study with the collection of both nasopharyngeal swab specimen and saliva specimen. The Antigen Tests (with Saliva specimens) and RT-PCR Tests (with Nasopharyngeal swab specimens) are tested in parallel. 134 samples of the symptomatic patients (within 2 days of onset) and 110 samples of healthy people were tested with Antigen Test and RT-PCR Test. The performance of the assay was confirmed with a series of clinical specimens by testing of specimens in a blinded fashion. The high sensitivity RT-PCR Test which uses a chemical lysis step followed by solid phase extraction of nucleic acid (e.g., silica bead extraction) was used as the comparator method.

Saliva		RT-PCR Test		Total
		Positive	Negative	
Antigen Test	Positive	123	0	123
	Negative	11	110	121
Total		134	110	244

Analysis	Statistics	95% Confidence Interval Score with Continuity Correction
Clinical Sensitivity	123 / 134 = 91.79%	(85.45% ~ 95.63%)
Clinical Specificity	110 / 110 = 100.00%	(95.79% ~ 99.92%)
True Positive (Positive Predict Value)	123 / 123 = 100.00%	(96.23% ~ 99.93%)
False Positive	0 / 123 = 0.00%	(0.00% ~ 3.77%)
True Negative (Negative Predict Value)	110 / 121 = 90.91%	(83.96% ~ 95.15%)
False Negative	11 / 121 = 9.09%	(4.85% ~ 16.04%)
Total Agreement Rate	(123 + 110) / 244 = 95.49%	(91.85% ~ 97.61%)

### LIMIT OF DETECTION

The final LOD is P5 = 1.25x10<sup>2</sup> TCID<sub>50</sub>/mL

No. of Reference	Concentration	Test
P1	2.0x10 <sup>5</sup> TCID <sub>50</sub> /mL	+ (20/20)
P2	2.0x10 <sup>4</sup> TCID <sub>50</sub> /mL	+ (20/20)
P3	0.5x10 <sup>3</sup> TCID <sub>50</sub> /mL	+ (20/20)
P4	0.25x10 <sup>3</sup> TCID <sub>50</sub> /mL	+ (20/20)
P5	1.25x10 <sup>2</sup> TCID <sub>50</sub> /mL	+ (20/20)
P6	0.625x10 <sup>2</sup> TCID <sub>50</sub> /mL	+ (17/20)
P7	0.3125x10 <sup>2</sup> TCID <sub>50</sub> /mL	+ (8/20)
P8	0.156x10 <sup>2</sup> TCID <sub>50</sub> /mL	+ (1/20)


### INTERFERENCE

Substance	Concentration	Negative	Positive	Interference (Yes / No)
Mucin	20 mg/mL	- (30/30)	+ (30/30)	No
Bilirubin	20 ug/mL	- (30/30)	+ (30/30)	No
Hemoglobin	15 mg/mL	- (30/30)	+ (30/30)	No
Rheumatoid Factor	20 mg/mL	- (30/30)	+ (30/30)	No
Antinuclear Antibodies	20 mg/mL	- (30/30)	+ (30/30)	No
Lipoproteins	20 mg/mL	- (30/30)	+ (30/30)	No
Glucose	20 mg/mL	- (30/30)	+ (30/30)	No
Uric Acid	200 ug/mL	- (30/30)	+ (30/30)	No
Acetaminophen	200 ug/mL	- (30/30)	+ (30/30)	No
Acetoacetic Acid	200 ug/mL	- (30/30)	+ (30/30)	No
Acetylsalicylic Acid	200 ug/mL	- (30/30)	+ (30/30)	No
Benzoyllecgonine	100 ug/mL	- (30/30)	+ (30/30)	No
Caffeine	200 ug/mL	- (30/30)	+ (30/30)	No
EDTA	800 ug/mL	- (30/30)	+ (30/30)	No
Ethanol	1.0%	- (30/30)	+ (30/30)	No
Gentisic Acid	200 ug/mL	- (30/30)	+ (30/30)	No
P- Hydroxybutyrate	200,000 ug/mL	- (30/30)	+ (30/30)	No
Methanol	10.0%	- (30/30)	+ (30/30)	No
Phenothiazine	200 ug/mL	- (30/30)	+ (30/30)	No
Phenylpropanolamine	200 ug/mL	- (30/30)	+ (30/30)	No
Salicylic Acid	200 ug/mL	- (30/30)	+ (30/30)	No

### CROSS-REACTIVITY

Pathogens/ Organisms	Negative	Positive	Cross-Reactivity (Yes / No)
MERS-coronavirus	- (30/30)	+ (30/30)	No
Human coronavirus 229E	- (30/30)	+ (30/30)	No
Human coronavirus OC43	- (30/30)	+ (30/30)	No
Human coronavirus NL63	- (30/30)	+ (30/30)	No
Human coronavirus HKU1	- (30/30)	+ (30/30)	No
Adenovirus	- (30/30)	+ (30/30)	No
Human Metapneumovirus	- (30/30)	+ (30/30)	No
Parainfluenza virus 1-4	- (30/30)	+ (30/30)	No
Influenza A & B	- (30/30)	+ (30/30)	No
Enterovirus	- (30/30)	+ (30/30)	No
Respiratory syncytial virus	- (30/30)	+ (30/30)	No
Rhinovirus	- (30/30)	+ (30/30)	No
Haemophilus influenzae	- (30/30)	+ (30/30)	No
Streptococcus pneumoniae	- (30/30)	+ (30/30)	No
Bordetella pertussis	- (30/30)	+ (30/30)	No
Mycoplasma pneumoniae	- (30/30)	+ (30/30)	No
Chlamydia pneumoniae	- (30/30)	+ (30/30)	No
Legionella pneumophila	- (30/30)	+ (30/30)	No
Mycobacterium tuberculosis	- (30/30)	+ (30/30)	No
Candida albicans	- (30/30)	+ (30/30)	No

Symbol			
	Reference		In Vitro Diagnosis
	Batch Code		Read the instructions before use
	Contains enough for <N> tests		Do not reuse Use only once
	Reserve : 2-30°C (35.6-86 °F)		Maker
	Date of Expiry		Danger
	European authority certification		CE Symbol

Lotus NL B.V.  
Koningin Julianaplein 10  
2595 AA's-Gravenhage

Careshield GmbH  
Veendyk 12c 46519 Alpen  
Registergericht: Amtsgericht Kleve  
Registernummer: HRB 16521

Ningbo Beautiful Life Medical Biotechnology Development Co., Ltd.  
Address: Room 305, Building 4, NO.7 Chuangye Avenue,  
Free Trade West Zone, Ningbo, China